

Every product arrives with a surface. Choose your preference:

- Stainless Steel 304/316: Durable and corrosion-resistant metal.
- Ceramic: Suitable for high temperatures and highly corrosive environments.
- Epoxy: Effective for use with solvents.
- HPL (High-Pressure Laminate): Resistant to solvents, petrol, and suitable for clean rooms.
- PP (Polypropylene): Resistant to corrosive substances.
- PPS (Polyphenylene Sulfide): Fire-retardant material, effective against corrosives.

Work Surface Material	Available colors	Surface thickness
Stainless Steel	Brushed silver, mirrored silver	1 mm (0.03 inch), 1.2 mm (0.04 inch), 1.5 mm (0.05 inch), 2 mm (0.07 inch)
Epoxy	Gray, white, blue, black	19 mm (0.7 inch)
HPL (High Pressure Laminate)	Gray, black	13 mm (0.5 inch), 16 mm (0.6 inch), 19 mm (0.7 inch)
Polypropylene	Gray, blue, white, black, ivory	8 mm (0.3 inch), 12 mm (0.4 inch), 15 mm (0.5 inch), 20 mm (0.7 inch)
PPS (Polyphenylene Sulfide)	Ivory	8 mm (0.3 inch), 12 mm (0.4 inch)

