



TOPAIR
SYSTEMS

Precision Stainless Steel Laboratory Solutions

TOPAIR





TopAir Systems

Engineered for Precision. Built for Performance.

About the Company

TopAir Systems is a leading company specializing in the design, development, and manufacturing of advanced solutions for controlled environments, including chemical laboratories, pharmaceutical and biotechnology industries, electronics, and cleanrooms.

The company combines advanced engineering expertise with years of industry experience to deliver high-quality products tailored to customer requirements, international standards, and evolving operational needs.

Stainless Steel Product Line

Our stainless steel product line is designed for demanding environments and includes a wide range of professional laboratory furniture and equipment.

All products are manufactured from high-quality stainless steel: **304 | 316 (customer selection)**

- High resistance to corrosion and chemicals
- Suitable for wet environments and cleanrooms
- Smooth finish for easy cleaning and sterilization
- Design that minimizes particle and contaminant accumulation

Customization & Engineering Flexibility

While this catalog presents standard product sizes, all products can be fully customized to meet specific requirements:

- Custom dimensions based on customer needs
- Structural and functional modifications
- Integration of accessories and add-ons
- Design and manufacturing of tailor-made solutions

Our engineering team supports each project from concept to delivery, ensuring precise and effective solutions for unique applications.

Advanced Material Integration

In addition to stainless steel products, the company also manufactures solutions made from: **Polypropylene (PP)**

- Combination of stainless steel and polypropylene in a single product
- Optimized solutions for chemical and functional requirements

TopAir Systems operates according to strict international standards and is certified with: **ISO 9001 | ISO 14001**

- Controlled manufacturing processes
- High product quality
- Commitment to environmental responsibility

Our Commitment

We are committed to delivering reliable, high-quality, and fully customized solutions, with a strong focus on innovation, precision, and professional service throughout every stage of the project.



For a customized product contact us: sales@topairsystems.com



Storage Cabinet

Functional for multiple cleanroom storage solutions



Name of Product	Model	Size (mm)		Size (in)	
		External	Internal	External	Internal
Single Door Storage Cabinet	S-1DO-800	800 x 500 x 700	795 x 475 x 567	31.5 x 19.69 x 27.56	31.3 x 18.7 x 22.32
Double Door Storage Cabinet	S-2DO-1200	1200 x 500 x 700	1195 x 475 x 567	47.2 x 19.69 x 27.56	47 x 18.7 x 22.32

For a customized product contact us: sales@topairsystems.com



Storage Cabinet with Drawers

Multifunctional cabinets designed for storage of hardware, medical devices, electronics and pharmaceuticals



Name of Product	Model	Size		Internal Dimensions	
		mm	in	mm	in
Storage Cabinet with Drawers	S-4DR- 600	600 x 600 x 800	23.62 x 23.62 x 31.5	450 x 500 x 150	17.71 x 19.69 x 5.90

For a customized product contact us: sales@topairsystems.com



Stainless Steel Cage

Used in dry environments, provides both strength and durability



Name of Product		Stainless Steel Cage
Model		S-CG- 1100
Size	mm	1100 x 1000 x 2175
	in	43.31 x 39.37 x 85.63
Shelf Size	mm	450 x 500 x 150
	in	17.71 x 19.68 x 5.90
Shelf Depth	mm	580
	in	22.83

For a customized product contact us: sales@topairsystems.com



Shoe Cabinet

Easy to-clean cabinet provides an organized area for clean room booties and shoes



Name of Product	Model	Size		Each space size	
		mm	in	mm	in
Shoe Cabinet	S-SC- 1000	1000 x 350 x 500	39.37 x 13.78 x 19.69	300 x 300 x 200	11.81 x 11.81 x 7.87

For a customized product contact us: sales@topairsystems.com



Lockers

Sloping top minimizes particle retention, offers minimal air disruption in gowning rooms with ceiling filter/blowers



Name of Product		Locker
Model		S-CLC 900
Size	mm	900 x 450 x 1950
	in	35.43 x 17.72 x 76.77
Each space size	mm	260 x 400 x 840
	in	10.24 x 15.75 x 33.07
Door size	mm	260 x 840
	in	10.24 x 33.07

For a customized product contact us: sales@topairsystems.com



Laboratory Cart

Provides basic storage and transport for dry environments



Name of Product	Model	Size		Shelf Dimensions	
		mm	in	mm	in
Laboratory Cart	S-CR-930	930 x 630 x 1000	36.61 x 24.8 x 39.37	930 x 630	36.61 x 24.8

For a customized product contact us: sales@topairsystems.com

Mesh Cart

Used for cleanroom materials



Name of Product	Model	Left		Righ		Shelf Dimensions	
		mm	in	mm	in	mm	in
Mesh Cart	S-WS-1200	600 x 400 x 750	23.62 x 15.75 x 29.53	600 x 400 x 400	23.62 x 15.75 x 15.75	600 x 400	23.62 x 15.75

For a customized product contact us: sales@topairsystems.com



Wire Shelving Shelf

Protect materials from sliding across shelves




Name of Product	Model	Size		Shelf Dimensions	
		mm	in	mm	in
Wire Shelving Shelf	S-WS-1200	1200 x 450 x 1600	47.24 x 17.72 x 62.99	1200 x 350	47.24 x 13.78

For a customized product contact us: sales@topairsystems.com



Work Table

Solid top, smooth surface, easy to clean and sterilize



Name of Product	Model	Size		Shelf Dimensions	
		mm	in	mm	in
Work Table	S-TA-1000	1000 x 500 x 800	39.37 x 19.69 x 31.5	1000 x 500	39.37 x 19.69

For a customized product contact us: sales@topairsystems.com

Table & Stool Set

Clean room table and stools set



Name of Product		Table & Stool Set
Model		TSS-001
Table size	mm	1000 x 600 x 800
	in	39.37 x 23.62 x 31.5
Stool size	mm	350 x 260 x 560
	in	13.78 x 10.24 x 22.05
Table shelf Dimensions	mm	1000 x 600
	in	39.37 x 23.62
Stool shelf Dimensions	mm	350 x 260
	in	13.78 x 10.24

For a customized product contact us: sales@topairsystems.com



Sink Station

Smooth surface, prevents accumulation of germs and other contaminants



Name of Product	Model	Size (mm)		Size (in)	
		External	Internal	External	Internal
Sink Station	S-STS-550	550 x 450 x 950	500 x 400 x 900	25.65 x 17.23 x 37.4	19.69 x 15.75 x 35.43

For a customized product contact us: sales@topairsystems.com



Clean Room Rack

Ideal for applications requiring solid shelving and maximum corrosion resistance




Name of Product		Clean Room Rack
Model		S-4CR-1200
Size	mm	1200 x 600 x 1800
	in	47.24 x 23.62 x 70.87
Shelf Dimensions	mm	1200 x 600
	in	47.24 x 23.62
Shelf Depth	mm	530
	in	20.87

For a customized product contact us: sales@topairsystems.com

Laboratory Ladder


Easy to reach, easy to clean step ladder



Name of Product	Model	Size		Shelf Dimensions	
		mm	in	mm	in
Laboratory Ladder	S-CLL -670	670 x 600 x 1545	26.38 x 23.62 x 60.83	670 x 600	26.38 x 23.62

For a customized product contact us: sales@topairsystems.com

Clean Bench






TOPAIR
SYSTEMS

More Solutions for Your Laboratory

Polypropylene fume hood integrated with a stainless steel base cabinet



For a customized product contact us: sales@topairsystems.com

TOPAIR



TOPAIR
SYSTEMS

www.topairsystems.com

Tel : 1-855-6-TOPAIR / International: +1-855-686-7247

Headquarters: 300 First Avenue, Suite 102, Needham, MA 02494 USA

Email: sales@topairsystems.com

